



PEACOCK LABORATORIES, INC. • 1901 S 54th St. • Philadelphia, PA 19143  
Telephone: 215-729-4400 • [www.peacocklabs.com](http://www.peacocklabs.com)

## **Technical Data**

### **Coppering Solutions**

#### **USES:**

Used for plating copper as a protective coating on top of the silver during the manufacture of mirrors by conveyORIZED mirror lines.

#### **Solution Chemistry:**

Consists of Peacock Copper Solution and Iron reducing powder. When diluted and used according our instructions, 20-22 mg of Copper is deposited per Sq. Ft of the silver mirror.

Adhesion: (measured by Tape Pull Test on silver + copper) - Nearly zero failure

#### **Packaging:**

Sold as one unit, consisting of one 5 gallon container of concentrated copper solution and 5 lbs. of Iron powder. Note: 1 unit will make 60 gallons of ready to spray solution.

#### **Price:**

Volume discounts are available. Call for pricing.

## **DISCLAIMER**

The above information was obtained from sources which we believe are reliable and correct, but does not purport to be inclusive and shall be used only as a guide. The information is provided without any representation or warranty expressed or implied regarding accuracy or correctness. The conditions or methods of handling, storage, use and disposal of the product are beyond our control and may be beyond our knowledge. For this and other reasons, Peacock Laboratories, Inc. and its employees do not assume responsibility and expressly disclaim liability for loss, damage, or expense arising out of, or in any way connected with, the handling, storage, use or disposal of the product.

**FOR INDUSTRIAL USE ONLY**